


TriMed Wrist Fixation System




# Bridge Plate









### **Exposure**

- Perform closed reduction with traction.
- Make an incision over the midshaft of either the 2<sup>nd</sup> or 3<sup>rd</sup> metacarpal.
- A second incision is made just proximal to the thumb outcropping muscles. Expose the radial shaft.
- An optional 3<sup>rd</sup> incision can be made over Lister's to mobilize the EPL or to create a graft portal.


## **Plate Insertion**

- Using the Bridge Plate Inserter, insert plate from proximal to distal, or from distal to proximal.
- Pass the plate under the extensor tendons along the surface of the bone until it appears through the opposite incision.

### **Provisional Distal Fixation**


- Position the plate over the desired metacarpal to allow placement of distal 2.7mm screws.
- Place a non-locking 2.7mm screw using a 2.0mm (yellow) drill in the most distal hole to lag the plate to bone.

# Bridge Plate



## **Secure Plate and Restore Length**

- With the wrist in neutral rotation, center the plate over the metacarpal and the radial shaft.
- Place a non-locking 3.2mm screw using a 2.3mm (red) drill into the distal end of the slotted hole.
- To gain additional length, loosen screw 1/4 turn, apply traction and retighten screw.



## **Final Fixation**

- Secure the plate by placing additional screws
  - 2.7mm locking or non-locking screws distally
  - 3.2mm locking or non-locking screws proximally


**Note:** When using locking drill guides or quick guides ensure installment and placement is concentric to the screw hole. Off-axis guide placement can result in screws not locking into the plate; locking screws can only be used on-axis.

#### TIPS

- Displaced or unstable volar ulnar corner, dorsal ulnar corner, or die-punch fragments may require augmented fixation.
- Buttress Pins<sup>™</sup>, Hook Plates<sup>™</sup>, or Pin Plates<sup>™</sup> can be used as adjunct fixation to a Bridge Plate.




All implants made from surgical grade stainless steel



#### **Bridge Plate**

BRGP



TriMed, Inc. / 27533 Avenue Hopkins / Santa Clarita, CA 91355 USA / 800-633-7221 / www.trimedortho.com

For indications, contraindications, warnings and precautions related to TriMed Wrist Fixation System reference IFU on trimedortho.com/ifu.

The technique presented is one suggested surgical technique. The decision to use a specific implant and the surgical technique must be based on sound medical judgment by the surgeon that

This document is controlled by TriMed, Inc. When downloaded, printed, and/or copied, this document becomes uncontrolled, and users should always check trimedortho.com for the latest version.

See trimedortho.com/patents for all patent information.

takes into consideration factors such as the circumstances and configuration of the injury.