

# Diaphyseal humerus

**135.1005**

## Diaphyseal humerus plate short

Length: 120 mm

Diaphysis holes	Locked: 8
	Free: 0




**135.1006**

## Diaphyseal humerus plate long

Length: 180 mm

Diaphysis holes	Locked: 10
	Free: 0



Scale 1:1

■ Autlocking screws  $\varnothing$  3.5 from 10 to 70 mm

### Materials

Plates: Stainless steel AISI 316 LVM - ISO 5832-1

Bushings and screws: Titanium Ti6Al4V - ISO 5832-3

# Diaphyseal tibia (and humerus)

**151.4100**

## ACP straight plate short

Length: 122 mm

**Diaphysis holes** Locked: 6  
Free: 2

**151.4101**

## ACP straight plate medium

Length: 169 mm

**Diaphysis holes** Locked: 6  
Free: 2

**151.4102**

## ACP straight plate long

Length: 209 mm


**Diaphysis holes** Locked: 8  
Free: 2

**151.4103**

## ACP straight plate extra-long

Length: 249 mm

**Diaphysis holes** Locked: 10  
Free: 2



**150.4500**

## Eyelet

To be used with Ø 2.0 mm cerclage wire



**150.4501**

## Threaded eyelet

To be used with Ø 2.0 mm cerclage wire

Scale 1:1.5



■ Autlocking screws Ø 5.0 from 12 to 110 mm

### Materials

Plates: Stainless steel AISI 316 LVM - ISO 5832-1  
Bushings and screws: Titanium Ti6Al4V - ISO 5832-3  
Eyelet: Titanium Grade 5

# Serie Mini/3.5

S30A - MINI/3.5 Instrument Case - 1 pcs



S21 - Ø 2.2 mm Conical Drill Guide - 1 pcs



S210 - Ø 2.2 mm Universal Drill Guide - 1 pcs



S213 - Ø 1.2 mm K-Wire Sleeve - 1 pcs



S215 - Ø 2.2 mm Conical Graduated Drill Guide - 1 pcs



S22 - Ø 2.2 Drill Bit, 160 mm - 2 pcs

S235 - K-Wire Ø 1.2 mm x 120 mm - 2 pcs



S24 - 2 mm Hex Screwdriver - 1 pcs



S246 - Compression Device for Locking Screws - 1 pcs



S25 - Bushing Carrier/Extractor - 1 pcs



S33 - Depth Meter - 1 pcs

S334 - K-Wire Ø 1.5 mm x 200 mm - 2 pcs

S335 - K-Wire Ø 1.5 mm x 120 mm - 2 pcs



S340 - 2.5 mm Hex Screwdriver - 1 pcs



S345 - Quick Coupling 2.5 mm Hex Shaft - 1 pcs



S346 - Compression Device for Locking Screws - 1 pcs



S35 - Bushing Carrier/Extractor - 1 pcs



S38 - Ø 2.8 mm Conical Drill Guide - 1 pcs



S382 - Ø 1.5 mm K-Wire Sleeve - 1 pcs



S384 - Ø 2.8 mm Conical Graduated Drill Guide - 1 pcs



S385 - Ø 2.8 mm Drill Guide for Cortical Screws - 1 pcs



S389 - Ø 2.8 mm Long Conical Drill Guide - 1 pcs



S39 - Ø 2.8 Drill Bit, 200 mm - 1 pcs



S390 - Ø 2.8 Drill Bit, 150 mm - 1 pcs



S394 - Ø 2.8 Drill Bit, 130 mm - 1 pcs



SVM - Ø 2.5 mm Temporary Stabilizing Screw - 2 pcs



SVS - Ø 3.5 mm Temporary Stabilizing Screw - 2 pcs

SET20/30 - Complete instruments set serie MINI/3.5