

# Excellence

## Proximal tibia

151.2101 (L) - 151.2102 (R)

### Proximal tibia plate short

Length: 83 mm

**Epiphysis holes** Locked: 4  
Free: 0

**Diaphysis holes** Locked: 2  
Free: 0

151.2103 (L) - 151.2104 (R)

### Proximal tibia plate medium

Length: 143 mm

**Epiphysis holes** Locked: 4  
Free: 0

**Diaphysis holes** Locked: 5  
Free: 0

151.2105 (L) - 151.2106 (R)

### Proximal tibia plate long

Length: 203 mm

**Epiphysis holes** Locked: 4  
Free: 0

**Diaphysis holes** Locked: 8  
Free: 0

151.2107 (L) - 151.2108 (R)

### Proximal tibia plate extra-long

Length: 263 mm

**Epiphysis holes** Locked: 4  
Free: 0

**Diaphysis holes** Locked: 11  
Free: 0



Scale 1:1.5

■ Autlocking screws  $\varnothing$  5.0 from 12 to 110 mm

#### Materials

Plates: Stainless steel AISI 316 LVM - ISO 5832-1

Bushings and screws: Titanium Ti6Al4V - ISO 5832-3

# Medial proximal tibia

151.2011

## Proximal tibia medial plate

Length: 124 mm

Epiphysis holes Locked: 3  
Free: 0

Diaphysis holes Locked: 4  
Free: 0



■ Autlocking screws  $\varnothing$  5.0 from 12 to 110 mm

### Materials

Plates: Stainless steel AISI 316 LVM - ISO 5832-1

Bushings and screws: Titanium Ti6Al4V - ISO 5832-3

Scale 1:1.5

# Serie 5

## Instruments Serie 5

S500 - 5/PFF - Instrument Tray - 1 pcs



S51 - Ø 4.0 mm Conical Drill Guide - 2 pcs



S512 - Ø 1.5 mm K-Wire Sleeve - 1 pcs



S514 - Ø 4.0 mm Conical Graduated Drill Guide - 1 pcs



S52 - Ø 4.0 mm Drill Bit, 200 mm - 1 pcs



S53 - Depth Meter - 1 pcs

S535 - K-Wire Ø 1.5 mm x 200 mm - 2 pcs



S540 - 3.0 mm Hex Screwdriver - 1 pcs



S55 - Bushing Carrier/Extractor - 1 pcs



S56 - Quick Coupling 3.0 mm Hex Shaft - 1 pcs



S562 - Screw/Bushing Disassembler - 1 pcs

## Dedicated instruments PFF

S102 - Threaded Guide Wire Ø 2.5 mm x 270 mm - 2 pcs

S59 - Ø 5.0 mm Conical Drill Guide - 2 pcs

S591 - 3.5 mm Hex Screwdriver - 1 pcs

S592 - PFF Depth Meter - 1 pcs

S593 - GT Hook Holder - 1 pcs

S594 - Ø 2.5 mm K-Wire Guide - 2 pcs

S595 - Ø 5.0 mm Cannulated Drill Bit, 200 mm - 1 pcs

SVL - Ø 5.0 mm Temporary Stabilizing Screw - 2 pcs

SET50/PFF - Complete instruments set serie 5 / PFF