

# PHF

## Proximal humerus


135.1203 (L) - 135.1204 (R)

### PHF Proximal humerus plate medium

Length: 82 mm

**Epiphysis holes** Locked: 7  
Free: 0

**Diaphysis holes** Locked: 2  
Free: 1



135.1205 (L) - 135.1206 (R)

### PHF Proximal humerus plate long

Length: 125 mm

**Epiphysis holes** Locked: 7  
Free: 0

**Diaphysis holes** Locked: 4  
Free: 2

135.1207 (L) - 135.1208 (R)

### PHF Proximal humerus plate extra-long

Length: 180 mm

**Epiphysis holes** Locked: 7  
Free: 0

**Diaphysis holes** Locked: 7  
Free: 2

■ Autlocking screws Ø 3.5 from 10 to 70 mm

■ Cortical screws Ø 3.5 from 10 to 70 mm

#### Materials

Plates and cortical screws: Stainless steel AISI 316 LVM - ISO 5832-1

Bushings and autlocking screws: Titanium Ti6Al4V - ISO 5832-3

Scale 1:1

# Serie Mini/3.5

S30A - MINI/3.5 Instrument Case - 1 pcs



S21 - Ø 2.2 mm Conical Drill Guide - 1 pcs



S210 - Ø 2.2 mm Universal Drill Guide - 1 pcs



S213 - Ø 1.2 mm K-Wire Sleeve - 1 pcs



S215 - Ø 2.2 mm Conical Graduated Drill Guide - 1 pcs



S22 - Ø 2.2 Drill Bit, 160 mm - 2 pcs

S235 - K-Wire Ø 1.2 mm x 120 mm - 2 pcs



S24 - 2 mm Hex Screwdriver - 1 pcs



S246 - Compression Device for Locking Screws - 1 pcs



S25 - Bushing Carrier/Extractor - 1 pcs



S33 - Depth Meter - 1 pcs

S334 - K-Wire Ø 1.5 mm x 200 mm - 2 pcs

S335 - K-Wire Ø 1.5 mm x 120 mm - 2 pcs



S340 - 2.5 mm Hex Screwdriver - 1 pcs



S345 - Quick Coupling 2.5 mm Hex Shaft - 1 pcs



S346 - Compression Device for Locking Screws - 1 pcs



S35 - Bushing Carrier/Extractor - 1 pcs



S38 - Ø 2.8 mm Conical Drill Guide - 1 pcs



S382 - Ø 1.5 mm K-Wire Sleeve - 1 pcs



S384 - Ø 2.8 mm Conical Graduated Drill Guide - 1 pcs



S385 - Ø 2.8 mm Drill Guide for Cortical Screws - 1 pcs



S389 - Ø 2.8 mm Long Conical Drill Guide - 1 pcs



S39 - Ø 2.8 Drill Bit, 200 mm - 1 pcs



S390 - Ø 2.8 Drill Bit, 150 mm - 1 pcs



S394 - Ø 2.8 Drill Bit, 130 mm - 1 pcs



SVM - Ø 2.5 mm Temporary Stabilizing Screw - 2 pcs



SVS - Ø 3.5 mm Temporary Stabilizing Screw - 2 pcs

SET20/30 - Complete instruments set serie MINI/3.5